



Application Instructions


Product: COOLFORCE Shoulder Brace
Model: 5510

- A quick inflation system that is easy to use, which can improve outcomes by managing pain and swelling.
- Removable ice pack.
- For use in the healing of sprains, sports injuries, arthritis and muscle fatigue.




Order# 5510 **Size** Universal

Sugg. HCPC: L3670



Picture 2



Picture 3



1. Check packaging first for (1) compressor (2) gel pack (3) one strap (4) pump with tube.
2. Refrigerate gel pack at least 2 hours before usage.
3. Put gel pack on the compressor (Picture 1) after refrigeration.
4. Connect the tube to the ON/OFF switch and turn switch to ON.
5. Pull the strap with black velcro through the square hook (from "IN" to "OUT") then fasten it (Picture 1).
6. Place the shoulder pad over your shoulder (A) and attach the velcro first Red then Blue on the proper position. (Picture 2). Picture 2 applies for the left shoulder. If for the right shoulder, use the opposite direction.
7. Pass the strap through your left or right armpit and pull it through the square hook (from "IN" to "OUT"). Then adjust velcro until comfortable (Picture 3)
8. If you wish to switch shoulders, loosen the velcro strap in front of chest. Remove the device and switch shoulders, then adjust the strap in front of chest until comfortable.
9. When turning ON/OFF switch, please hold it with one hand and turn it with the other hand to ensure the switch is on the appropriate position, then pump.