

At Advanced Orthopaedics, a new day of healing has arrived. Our orthopedic braces are scientifically engineered to accelerate the healing process and to reduce the body's recovery time. Each product is carefully designed to protect the body while providing superior comfort and support.





Welcome to Advanced Orthopaedics®

#### **Customer Service**

Our friendly and knowledgeable Customer Service Specialists are available from 8:00 am - 5:00 pm EST to take your order or answer any questions about our products.

#### Ordering

Phone 800.270.7074 215.256.9006 Fax

Email orders@advanced-orthopaedics.com

Online advanced-orthopaedics.com

#### Shipping

We offer same day shipping. All orders are shipped FedEx & USPS ground unless specified otherwise. Orders are shipped directly from our warehouse located in Bloomsburg, PA.

#### Warranty

We guarantee complete satisfaction on all of our products. Only the finest grade materials are used to assure superior quality and durability.

If for any reason you are unsatisfied with our products' quality and workmanship, please contact our Customer Care department.





## **PRODUCT CATEGORIES**

Cervical, Spine & Abdominal	4
Knee	13
Wrist & Elbow	23
Walking Boots & Night Splints	32
Ankle & Foot	37
Shoulder	46
Therapeutic Aides	50
Fracture Management	53
Compression Hosiery	58

This icon indicates that the product is available with Retail Packaging.




These QR codes can be scanned by your smart phone and directs you to the Product Detail Page.









#### **Universal Cervical Collar**

- Styled to give uniform support
- · Medium-firm density foam covered with tubular knit fabric
- Velcro <sup>™</sup> closure

**Indications:** Designed to help with whiplash, herniated cervical discs and other minor neck injuries. | Suggested HCPC L0120

Order #	Size	Length	Height
2000	Universal	18"	2½"
2030	Universal	18½"	3" Services
2035	Universal	21½"	3½"
2040	Universal	22"	4"



Height refers to measurement under the chin.



#### **Premium Cervical Collar**

- · Contoured to give uniform support
- Medium-firm density foam covered with tubular knit fabric
- Velcro ™ closure

Indications: Designed to help with whiplash, herniated cervical discs and other minor neck injuries. | Suggested HCPC L0120

Order #	Size	Length	Height
2003	Small	14½"	2½"
2005	Medium	16½"	3" <b>• ****</b>
2007	Large	18½"	3½"
2008	X Large	20½"	4" <b>10.000</b>





#### **EZ Glide Cervical Collar Series**

- 12" height adjustment settings for a custom fit on the front, ergonomic dynamic knob
- Adjustable b ack panel for a custom fit (EZ-180 model)
- · Larger tracheal opening

**Indications:** Designed to help with whiplash, herniated cervical discs and other minor neck injuries.

Order #	Size	Suggested	HCPC
EZ-172	Universal	L0172	<b>□%€</b> %□
EZ-174	Universal	L0174	
EZ-180	Universal	L0180	





#### **AO Cervical Collar**

- · Provides rigid immobilization, excellent for emergency use
- Reduces rotation and slightly hyperextends the cervical spine
- · Molded and shaped to conform to any patient
- Constructed of a lightweight closed cell foam material
- Velcro<sup>™</sup> closure facilitates easy application and removal
- Height refers to measurement under the chin

Indications: Moderate control of unwanted neck flexion /extension. Treat mid-neck strains and sprains. | Suggested HCPC L0172

Order #	Size	Height	Neck Circumference
02S	Small	21/4"	10" - 13"
02M	Medium	21/4"	13" - 16"
02L	Large	21/4"	16" - 19"
03S	Small	31/4"	10" - 13"
03M	Medium	31/4"	13" - 16"
03L	Large	31/4"	16" - 19"
04S	Small	41/4"	10" - 13"
04M	Medium	41/4"	13" - 16"
04L	Large	41/4"	16" - 19"
05S	Small	51/4"	10" - 13"
05M	Medium	51/4"	13" - 16"
05L	Large	51/4"	16" - 19"





#### Trac Collar™

- Portable and light weight traction collar gently stretches and relaxes muscles for pain relief
- · Fastens easily with buckles
- Includes carrying case
- Constructed of soft, plastic foam; internal inflatable bladders contain natural rubber latex

Indications: Used to reduce pain and discomfort of chronic and degenerative neck ailments. Not for recent sprains, strains, injuries or acute conditions. | Suggested HCPC E0856

Order #	Size	Neck Circumference
2050	Regular   Red	14" - 16"
2051	Large   Blue	16" - 18"





#### Transformer<sup>2</sup> Back Brace<sup>®</sup>

- Updated model with double pull design
- · Rigid panels provide superior posterior and anterior support
- Patient friendly, ensures greater patient compliance
- · Quick double strap pull controls your support level
- · Easy to put on and off
- · Breathable mesh belt
- · Produces intracavitary pressure to reduce load on invertebral discs

**Indications:** Spinal Stenosis, Spondylolisthesis, Spondylolysis, Chronic Back Pain, Spinal Instability, Post-Op Disk Surgery, Facet Syndrome, Degenerative, Hemiated, Bulging Disc, Post-Operative Lumber Fusion.

The easy to use "Slide-Out" adjustment system allows for a quick custom fit, making the Transformer<sup>2</sup> a truly universal back brace.

Suggested HCPC for 8027 is L0627 and L0642 Suggested HCPC for 8031 is L0631 and L0648 Suggested HCPC for 8037 is L0637 and L0650

The Transformer<sup>2</sup> has a universal belt size which fits Small, 26" waist circumference, through 4X Large, 60" waist circumference.



#### Transformer<sup>2</sup> Series

Order #	Size	Back Panel Height	Side panels
8027	Small - 4X Large	10"	No
8031	Small - 4X Large	14"	No
8037	Small - 4X Large	14"	Yes



#### Weave® Series

- · Single pull strap design
- Quick pull 6.1 mechanical advantage compression system
- Panels are sculpted to hug the contours of the back
- Rigid panels provide superior posterior and anterior support
- Light Weight
- Low profile sides conform naturally to different patients anatomies

Indications: Spinal Stenosis, Spondylolisthesis, Spondylolysis, Chronic Back Pain, Spinal Instability, Post-Op Disk Surgery, Facet Syndrome, Degenerative, Hemiated, Bulging Disc, Post-Operative Lumber Fusion.

Weave® 27Weave® 31Weave® 77Suggested HCPCSuggested HCPCSuggested HCPCL0627 and L0642L0631 and L0648L0637 and L065010" high back panel14" high back panel14" high back panel plus side panels for extra lateral stability

Weave® 27	Weave® 31	Weave® 77		
Order #	Order #	Order #	Size	Waist Circumference
733	743	773	Small	28" - 32"
735	745	775	Medium	32" - 36"
737	747	777	Large	36" - 40"
738	748	778	X Large	40" - 46"
739	749	779	2X Large	46" - 50"
739-9	749-9	779-9	3X Large	50" - 56"
739-99	749-99	779-99	4X Large	56" - 62"





#### APEX LSO<sup>2</sup> & APEX TLSO Back Brace®

- · Easy to use universal sizing system and adjustable handles
- · Adjustable pull handles to allow optimal compression
- · Breathable fabric allows air flow to keep patients cool and dry over long periods of time

Indications: Spinal Stenosis, Spondylolisthesis, Spondylolysis, Chronic Back Pain, Spinal Instability, Post-Op Disk Surgery, Facet Syndrome, Degenerative, Herniated, Bulging Discs, Post-Op Lumber Fusion.

Order #	Order #	Size	Suggested HCPC
APEX LSO <sup>2</sup> Back Brace 27/42	5027	Universal	L0627 and L0642
APEX LSO <sup>2</sup> Back Brace 31/48	5031	Universal	L0631 and L0648
APEX LSO <sup>2</sup> Back Brace 37/50	5037	Universal	L0637 and L0650
APEX LSO <sup>2</sup> Back Brace 39/51	5039	Universal	L0639 and L0651
APEX TLSO Back Brace 56/57	5056	Universal	L0456 and L0457



#### Ice Packs

- · Breathable, cool wicking material
- · Available in two sizes and attaches to the brace
- · Designed to fit a variety of braces in order to speed up recovery times and reduce pain

Order#	Size
27	Small (7.5 W x 6.5 T) Large (9.5 W x 13.5 T)
37	Large (9.5 W x 13.5 T)





### Lumbo Lite® Lumbar Orthosis

- · Streamline design allows for exact fit and comfort
- Posterior and anterior padded panels (rigid and moldable) which provide lower back and abdominal compression

Indications: Spinal Stenosis, Spondylolisthesis, Spondylolysis, Chronic Back Pain, Spinal Instability, Post-Op Disk Surgery, Facet Syndrome, Degenerative, Hemiated, Bulging Disc, Post-Operative Lumber Fusion. | Suggested HCPC L0627 and L0642

Order # Black	Order # White	Size	Waist Circumference
413 - B	413 - W	Small	26" - 32"
415 - B	415 - W	Medium	30" - 37"
417 - B	417 - W	Large	34" - 41"
418 - B	418 - W	X Large	38" - 47"
419 - B	419 - W	2X Large	44" - 54"
420 - B	420 - W	3X Large	50" - 62"
421 - B	N/A	4X Large	62" - 72"







### T.L.S.O. Support

- · Two flexible metal stays
- · Adjustable metal stays provide better support
- · Comfortable to wear and unnoticeable under clothing
- · Reinforcement straps provide extra support

Indications: Produces intercavitary pressure to reduce load in the intervertebral discs with rigid stays or panels. | Suggested HCPC L0450

Order #	Size	Waist Circumference
1093	Small	30" - 34"
1095	Medium	34" - 38"
1097	Large	38" - 42"
1098	X Large	42" - 46"
1099	2X Large	46" - 50"
1100	3X Large	50" - 54"



### EZ Pull Adjustable Back & Abdominal Support®

- · Reduces stress and strain on the lower back
- Lightweight, breathable material & universally-sized
- · One-Handed fully adjustable pulley fitting system
- · State-of-the-art sequential compression
- · Ideal for everyday use and can be worn under clothing

**Indications:** Universal size, designed to reduce stress and strain on the lower back. Lower back support, the lower back muscles are supported by the rigid back panels. Relieves back pain, weakness and fatigue. Easy single pull fitting system. | Suggested HCPC L0626

Order # 450	Size
450	Universal Size (28" - 52" Waist)







#### Low Profile SI Belt

- Stabilization of the sacroiliac (SI) joint through direct circumferential compression
- Mechanical advantage pulley system with 3 to 1 ratio
- · Customized stabilization pressure specific to patient

Indications: For post partum symphysis pubis dysfunction with pads removed. | Suggested HCPC L0621

Order #	Size	Waist Circumference
632	Small	Up to 35"
633	Standard	35" - 50"







### Premium Lumbar Sacral Support with Abdominal Belt

- Form-fitting 10" all-elastic with double pull side panels for additional support
- High quality Velcro<sup>™</sup> closure for durability
- Four flexible stays in back panel for excellent support

Indications: Ideal for providing support and compression for strains and sprains of the lower back. | Suggested HCPC L0625

Order #   White	Order #   Black	Size	Waist Circumference
501 - W	501 - B	X Small	20" - 24"
503 - W	503 - B	Small	24" - 30"
505 - W	505 - B	Medium	30" - 34"
507 - W	507 - B	Large	34" - 38"
508 - W	508 - B	X Large	38" - 42"
509 - W	509 - B	2X Large	42" - 46"
510 - W	510 - B	3X Large	46" - 50"
511 - W	511 - B	4X Large	50" - 54"
512 - W	512 - B	5X Large	54" - 58"





Larger sizes available upon request.

### Breathe Lite Back Support®

- Designed with a unique ventilated elastic construction making it cool and breathable
- Provides support and compression to the lower lumbar region and abdominal area
- Tapered front end period. Designed with an 8" high back panel, which is ideal for shorter or geriatric patients
- · Available in white and black

**Indications:** Ideal for providing support and compression for strains and sprains of the lower back.

Order #   White	Order #   Black	Size	Waist Circumference
573 - W	573 - B	Small	25" - 29"
575 - W	575 - B	Medium	28" - 32"
577 - W	577 - B	Large	31" - 35"
578 - W	578 - B	X Large	34" - 38"
579 - W	579 - B	2X Large	37" - 41"
580 - W	580 - B	3X Large	40" - 44"

Larger sizes available upon request.



### Premium Criss-Cross Lumbar Sacral Support

- Lightweight with classic criss-cross reinforcement design
- Elastic construction with four stays for added support
- 9" in back tapers to 4" in the front

Indications: Ideal for providing support and compression for strains and sprains of the lower back.

Order #	Size	Waist Circumference
513	Small	28" - 31
515	Medium	31" - 34"
517	Large	34" - 37"
518	X Large	37" - 40"
519	2X Large	40" - 43"
520	3X Large	43" - 46"





Larger sizes available upon request.

#### **Industrial Back Support**

- · Form fitting, lightweight material with double pull side panels for additional support
- Provides excellent lower back and abdominal support
- Suspenders adjust for a comfortable fit
- High quality Velcro<sup>™</sup> closure for durability

Indications: Ideal for providing support and compresion for stra sprains of the lower back.

Order #	Size	Waist Circumference	
403	Small	24" - 30"	
405	Medium	30" - 34"	
407	Large	34" - 38"	
408	X Large	38" - 42"	
409	2X Large	42" - 46"	
410	3X Large	46" - 50"	
411	4X Large	50" - 54"	■ 188

Larger sizes available upon request.





#### **Premium Abdominal Binder**

- · Available in 3" and 4" panels
- · Form-fitting, breathable abdominal binders conform to different patient shapes for improved compression
- · Made of soft foam and high quality elastic
- Velcro<sup>™</sup> closure
- · Latex free

Indications: Commonly used to augment the recovery process after abdominal surgery like exploratory laparotomy, cesarean section, bariatric surgery, hysterectomy or spinal surgery.

Order #	Size	Waist Circumference
2501	3 Panel (9" width)	30" - 45"
2502	3 Panel (9" width)	45" - 62"
2503	4 Panel (12" width)	30" - 45"
2504	4 Panel (12" width)	45" - 62"







### **Universal Rib Belt Support**

- · Constructed of 6" wide breathable elastic material
- · Female belt has contours around bust line and male belt is straight

Indications: Helps stabilize rib and sternum fractures by limiting chest expansion through compression.

Order #	Size	Rib Circumference
2400	Male	24" - 60"
2410	Female	24" - 60"







#### Premium Hernia Belt

- · Removable and adjustable pads
- Discreet
- · Double belt

Indications: Helps to reduce pain and discomfort of hernia bulging by applying resistive pressure over the location. | Suggested HCPC L8310

Order #	Size	Waist Circumference
493	Small	30" - 35"
495	Medium	35" - 41"
497	Large	41" - 46"











#### Deluxe Hinged Knee Brace

- · Medial and lateral hinge stays provide excellent support
- · Enclosed tubular donut buttress
- · Neoprene coated with 4-way stretch nylon fabric
- · Piping along thigh and calf areas of support to stop flexing
- · Tapered in back for comfortable fit
- · Easy removable hinges for custom fitting

Indications: Mild to moderate joint pain, sprain, strains of the knee. Provides support, stabilization and compression of the knee joint following ligament, tendon or meniscus injuries | | Suggested HCPC | 1810 and | 181

IIICIIISCUS	ilijulies. Į i	ouggested flor C L 1010	and Lion	_	
		Knee Circumference			
903	Small	12½" - 14½"	908	X Large	17½" - 20"

903	Small	12½" - 14½"	908	X Large	17½" - 20"		▣
905	Medium	14" - 16"	909	2X Large	19½" - 22" 21½" - 24"	-	Ę
907	Large	15½" - 18"	910	3X Large	21½" - 24"		
						0.00	

Measure around the center of the knee with leg extended.

### Deluxe Airprene Hinged Knee Brace

- · A unique perforated neoprene-like material that retains heat, but not water.
- · The perspiration generated by exercise flows through the AIRPRENE to keep the skin dry and comfortable.
- · Piping along thigh and calf areas of support to stop flexing
- · Tapered in back for comfortable fit
- · Easy removable hinges for custom fitting

Indications: Mild to moderate joint pain, sprain, strains of the knee.

Ligament o	r patellar in	istabilities.   Suggest	ed HCPC L1	810 and L1	812
Order #	Size	Knee Circumference	Order #	Size	Knee Circum

903 - AP	Small	12½" - 14½"
905 - AP	Medium	14" - 16"
907 - AP	Large	15½" - 18"

Measure around the center of the knee with leg extended.







### Wrap-Around Hinged Knee Brace

- · Excellent choice for patients who are not able to pull a knee sleeve brace over their feet
- · Medial and lateral hinge stays provide excellent support
- · Unique wrap-around design provides for an excellent custom fit
- · Elastic piping throughout the brace for better durability
- 12" length

Indications: Mild to moderate joint pain, sprain, strains of the knee. Provides support, stabilization and compression of the knee joint following ligament, tendon or meniscus injuries. | Suggested HCPC L1820

meniscus injunes.   Guggested Flor G E 1020						
Order #	Size	Knee Circumference	Order #	Size	Knee Circum	nference
603	Small	12½" - 14½"	609	2X Large	19½" - 22"	
605	Medium	14" - 16"	609-9	3X Large	21" - 24"	Man Brook
607	Large	15½" - 18"	609-99	4X Large	25" - 28"	
608	X Large	17½" - 20"	609-999	5X Large	29" - 32"	









#### Gator Wrap Universal Hinged Knee Brace®

- Excellent choice for patients who are not able to pull a knee sleeve brace over their feet
- · Medial and lateral hinge straps provide excellent support
- Unique wrap design provides for an excellent custom fit

Indications: Indicated for mild to moderate ligament joint pain, sprain, strains of the knee ligaments. Provides support, stabilization and compression of the knee joint following ligament or meniscus injuries. | Suggested HCPC L1820

Order # | Size 6130

Universal





## Min-Knee Hinged Knee Brace®

- · Wrap-around design allows for exact fit
- · Elastic piping throughout the brace for better durability
- · Medial and lateral hinge stays provide excellent support
- Neoprene coated with 4-way stretch nylon fabric
- · Tapered in back for comfortable fit
- · Easy removable hinges for custom fitting
- Polycentric full range of motion hinge with stops from 0° to 90°

Indications: Indicated for moderate to severe ligament joint pain, sprain, strains of the knee ligaments. Provides support, stabilization and compression of the knee joint following ligament or meniscus injuries. | Suggested HCPC L1832 and L1833

Order #	Size	Knee Circumference	
623	Small	12½" - 14½"	
625	Medium	14" - 16"	
627	Large	15½" - 18"	
628	X Large	17½" - 20"	
629	2X Large	19½" - 22"	
629-9	3X Large	21½" - 24"	
629-99	4X Large	25" - 28"	
629-999	5X Large	29" - 32"	

Measure around the center of the knee with leg extended.





### Cross-Fit Universal Hinged Knee Brace®

- · Medial and lateral hinge stays provide excellent support
- · Enclosed tubular donut buttress
- Neoprene coated with 4-way stretch nylon fabric
- · Piping along thigh and calf areas of support to stop flexing
- · Tapered in back for comfortable fit
- Polycentric full range of motion hinge with stops from 0° to 90°
- · Easy removable hinges for custom fitting

Indications: Indicated for moderate to severe ligament joint pain, sprain, strains of the knee ligaments. Provides support, stabilization and compression of the knee joint following ligament or meniscus injuries. | Suggested HCPC L1832 and L1833



Size Universal











### TM Wrap-Around Hinged Knee Brace®

- · The wrap-around design allows for exact fit and comfort
- · Provides mild to moderate medial and lateral support; in addition to the support offered by the strapping
- · Polycentric hinge provides extension
- Stops from 0° to 90°

Indications: Mild injuries of the meniscus mild instability, patellar retinaculum injuries, post op knee rehabilitation, mild sprains of the medial or lateral collateral ligaments. | Suggested HCPC L1832 and L1833

Order #	Size	Thigh Circumference
933	Small	Less than 18"
935	Medium	18" - 21"
937	Large	21" - 24"
938	X Large	24" - 27"
939	2X Large	27" - 30"
939-9	3X Large	30" - 35"



Measure around the center of the knee with leg extended

## Advanced Hinged Range of Motion Knee Brace

- · Telescoping thigh and calf for a custom fit
- · Easy to operate drop lock to ensure secure positioning
- · Polycentric hinge allows controlled range of motion with easy to adjust pin stops
- · Flexion / extension adjustments can be made easily
- · Universal brace is easy to fit

Indications: Mild injuries of the meniscus mild instability, patellar retinaculum injuries, post op knee rehabilitation, mild sprains of the medial or lateral collateral ligaments. | Suggested HCPC L1832 and L1833

Order#	Size	
917 - T	Universal	



#### **ACL Knee Brace**

- · Strong, lightweight composite frame.
- · User friendly ROM hinge for quick and easy adjustment

Indications: Provides stability and support following injury or reconstruction of the ACL, PCL, MCL, LCL and menisci. Can be used prophylactically to prevent hyperextension injuries. | Suggested HCPC L1845 and L1852

Order # Right	Order # Left		
843 - R	853 - L	Small	15½"
845 - R	855 - L	Medium	18
847 - R	857 - L	Large	21" - 23½"
848 - R	858 - L	X Large	23½" - 26½"
849 - R	859 - L	2X Large	2



Note: Thigh circumference should be measured approximately 6" above the knee.



#### OA Forte Unloader Knee Brace

- · Ultra-lightweight design and simple to apply
- · Easy to adjust straps with push-button buckles
- · Adjustable unloading hinge
- · Valgus and varus adjustment to fit individual's alignment

Indications: Ideal for patients with mild to moderate osteoarthritis of the knee. | Suggested HCPC L1843 and L1851

of the knee.   Sugge	ested HCPC L1843 and L1851	
Order#	Location	
960 - R 960 - L	Fits Right Lateral and Left Medial Fits Left Lateral and Right Medial	7000 E



### Panther Unloader Knee Brace®

- · Extremely light-weight
- Low-profile brace provides 3 points of leverage to reduce pressure inside the knee joint
- · Ideal for unilateral compartment OA knee users
- · Valgus and varus adjustment to fit individual's alignment

Indications: Ideal for patients with mild to moderate osteoarthritis of the knee. | Suggested HCPC L1843 and L1851

Order #	Location
950 - R	Fits Right Lateral and Left Medial
950 - 1	Fits Left Lateral and Right Medial



New

**Product** 



#### Fast-Fit Sleeve®

- Perfect addition for patients that wear hinged knee braces
- Lightweight cotton/lycra sleeve helps prevent braces from sliding and moving
- · Neoprene lining provides extra comfort to patient
- For use with hinged knee braces to minimize migration

Indications: For use with hinged knee braces. | Suggested HCPC L2397









### Sport Lite Knee Brace

- Enclosed neoprene tubular donut buttress
- · Two medial and two lateral spiral steel stays
- · Tapered in back for comfortable fit
- · Provides warmth and equal compression

Indications: Mild to moderate ligament joint pain, strains, sprains of the knee ligaments or patellar instabilities. | Suggested HCPC L1828

Order #	Size	Knee Circumfe	rence	
813	Small	12½" - 14½"		
815	Medium	14" - 16"		
817	Large	15½" - 18"		
818	X Large	17½" - 20"	5	
819	2X Large	19½" - 22"	3	
820	3X Large	21½" - 24"	No.	

Measure around the center of the knee with leg extended.



- "Perforated Neoprene" to allow moisture to escape skin area
- Two medial and two lateral spiral steel stays
- · Universal buttress pad for easy adjustment and exact fit
- · Tapered in back for comfortable fit
- · Piping along thigh and calf areas of support to stop flexing

**Indications:** Provides stability and support following injury or reconstruction of the ACL,PCL, MCL, LCL and menisci. Can be used prophylactically to prevent hyperextension injuries. | Suggested HCPC L1825

Order #	Size	Knee Circumference
803	Small	12½" - 14½"
805	Medium	14" - 16"
807	Large	15½" - 18"
808	X Large	17½" - 20"
809	2X Large	19½" - 22"

Measure around the center of the knee with leg extended.

#### Visco-Gel Silicone Knee Brace

- Flat knitted, four-way stretch elastic provides comfortable moderate compression
- Contoured silicone ring redistributes patella pressure
- Two medial and two lateral spiral steel stays

Indications: For tender knees resulting from minor sprains and strains. Suggested HCPC L1800  $\,$ 

Order #	Size	Knee Circumfe	rence
3333	Small	12" - 13½"	
3355	Medium	13½" - 15"	
3377	Large	15" - 16½"	
3388	X Large	16½" - 18"	
3399	2X Large	18" - 19½"	

Measure around the center of the knee with leg extended.







#### Universal Wrap Around Knee Brace

- Excellent choice for those patients who are not able to pull a knee sleeve over their feet
- Unique wrap-around design fits a wide range of sizes
- Two medial and two lateral spiral steel stays
- · Universal buttress pad for easy adjustment and exact fit

Indications: Mild to moderate ligament joint pain, strains, sprains of the knee ligaments or patellar instabilities. | Suggested HCPC L1800

Order #	Size	Knee Circumference
600	Universal	12" - 23"



Measure around the center of the knee with leg extended.

### Sport Knee Sleeve Support

- Neoprene coated with 4-way stretch nylon fabric
- · Tapered in back for comfortable fit
- · Open patella relieves pressure on the patella area

Indications: Mild to moderate ligament joint pain, strains, sprains of the knee ligaments or patellar instabilities. | Suggested HCPC L1825

Order #	Size	Knee Circumference	
303	Small	12½" - 14½"	
305	Medium	14" - 16"	
307	Large	15½" - 18"	
308	X Large	17½" - 20"	
309	2X Large	19½" - 22"	鳞
309-9	3X Large	21½" - 24"	僡

Measure around the center of the knee with leg extended.

#### Airprene Knee Sleeve

- Unique perforated neoprene-like material that retains heat, but not water
- The perspiration generated by exercise flows through the Airprene to keep the skin dry and comfortable
- Airprene helps the brace stay firmly in place by dissipating excessive moisture, eliminating the slipping and movement of the conventional neoprene brace

**Indications:** Mild to moderate ligament joint pain, strains, sprains of the knee ligaments or patellar instabilities. | Suggested HCPC L1825

Order #	Size	Knee Circumference	
303 - AP	Small	12½" - 14½"	
305 - AP	Medium	14" - 16"	
307 - AP	Large	15" - 18"	回线等数回
308 - AP	X Large	17½" - 20"	
309 - AP	2X Large	19½" - 22"	

Measure around the center of the knee with leg extended.









### Elastic Slip-On Knee Support (Open Patella)

- · Easy slip-on design
- · Designed to give measured support to a weak knee
- · Patella opening helps position sleeve properly
- 11" sleeve

Indications: Mild sprain or strain of the knee tendons or ligaments.

Order #	Size	Knee Circumfe	rence	
323	Small	12½" - 14½"		
325	Medium	14" - 16"	Say Super	回線機画
327	Large	15½" - 18"	1	
328	X Large	17½" - 20"	<b>8</b> ) Sec.	
3	323 325 327	Small Medium Large	Small 12½" - 14½"  Medium 14" - 16"  Large 15½" - 18"	Small 12½" - 14½"  Medium 14" - 16"  Large 15½" - 18"

Measure around the center of the knee with leg extended.



### Elastic Slip-On Knee Support (Closed Patella)

- · Easy slip-on design
- · Durable knitted elastic provides support with flexibility
- · Seamless construction will not irritate skin
- 11" sleeve

Indications: Mild sprain or strain of the knee tendons or ligaments.

Order #	Size	Knee Circumference		
333	Small	12½" - 14½"		
335	Medium	14" - 16"	District Control	
337	Large	15½" - 18"		
338	X Large	17½" - 20"	Ø 800	可能被保护

Measure around the center of the knee with leg extended.



#### Thermoskin Elastic Knee Stabilizer

- · Anatomically designed for comfort and fit
- Side stabilizers provide increased medial and lateral support to weak and injured knees
- · Available in retail packaging

Indications: Mild sprain or strain of the knee tendons or ligaments.

Order #	Size	Knee Circur	mference
R - 83646	Small	12" - 13½"	
R - 84646	Medium	14" - 15½"	
R - 85646	Large	16" - 17¾"	
R - 86646	X Large	18" - 19¾"	

Measure slightly underneath the knee cap.



#### Patella Knee Strap

- · Strap functions dynamically as knee flexes and extends
- · Easy to apply and adjust

Indications: Reduces pain for common joint conditions including Arthritis, Tendonitis and Chondromalacia. Ideal for minor sprains and strains of the knee ligaments.

Order #	Size
610	Universal







### Neoprene Thigh Sleeve Support

- · Designed for relief and protection against injured quadriceps and hamstring muscles
- · Provides warmth and equal compression
- · Excellent Support

Indications: Mild sprain, strain of quadriceps muscle.

Order #	Size	Thigh Circumference	
300	Small	18" - 20"	
300	Medium	20" - 22"	
300	Large	22" - 24"	
300	X Large	24" - 26"	

Measure around center of thigh.



### Neoprene Calf Sleeve Support

- · Designed for relief and protection against injured calf muscles
- · Provides warmth and equal compression
- · Excellent Support

Indications: Mild sprain, strain of gastrocnemius muscle.

Order #	Size	Calf Circumference	
301	Small	13" - 14"	
301	Medium	14" - 15"	回路線回
301	Large	15" - 17"	
301	X Large	17" - 19"	

Measure around the center of the calf.





#### Premium Sized Knee Immobilizer

- · Lightweight, cool and comfortable
- · Provides support and stabilization for the leg and knee
- · Three posterior "I" beam stays for extra rigidity
- · Rigid medial and lateral stays for greater immobilization
- Pro-pad insert included

Indications: Sprains, strains of the knee ACL, PCL, MCL, LCL injuries. | Suggested HCPC L1830

Order #	Size	Thigh Circumference
7003	Small	15" - 17"
7005	Medium	17" - 19"
7007	Large	19" - 21"
7008	X Large	21" - 23"
7009	2X Large	23" - 25"







Measure 4" above the patella.

#### Tri-Panel Knee Immobilizer

- · Tri-Panel design for more anatomical fit
- · Movable side panels for medial and lateral stay positioning
- Three contoured posterior stays
- · Two medial and two lateral stays
- Velcro<sup>™</sup> closure

Indications: Sprains, strains of the knee ACL, PCL, MCL, LCL injuries. | Suggested HCPC L1830

Order #	Size	Length	
703 - H	Universal	14"	
705 - H	Universal	20"	
707 - H	Universal	24"	

Each length is universal size.







#### **Premium Wrist Brace**

- · Ideal for carpal tunnel syndrome, wrist sprains, strains or after cast removal
- Maximum contouring at the palmar crease allows for full finger function
- Removable and adjustable palmar stay for maximum support
- Unique stockinette liner provides comfort and easy slip on application
- Velcro<sup>™</sup> closure for easy adjustment
- 8" length

Indications: Carpal tunnel syndrome, sprains, strains, tendonitis, wrist immobilization and post cast healing. | Suggested HCPC L3908

Order # Left	Order # Right	Size	Wrist Circumference
431 - L	421 - R	X Small	4½" - 5½"
433 - L	423 - R	Small	5½" - 6½"
435 - L	425 - R	Medium	6½" - 7½"
437 - L	427 - R	Large	7½" - 8½"
438 - L	428 - R	X Large	8½"+





### Lycra Lined Premium Wrist Brace

- · Soft interior for exceptional comfort
- · Removable palmer stay to allow correct fit
- · Lycra lining wicks moisture away from skin
- 8" length

Indications: Carpal tunnel syndrome, sprains, strains, tendonitis, wrist immobilization and post cast healing. | Suggested HCPC L3908

Order # Left	Order # Right	Size	Wrist Circumference
431 - L - C	421 - R - C	X Small	4½" - 5½"
433 - L - C	423 - R - C	Small	5½" - 6½"
435 - L - C	425 - R - C	Medium	6½" - 7½"
437 - L - C	427 - R - C	Large	7½" - 8½"
438 - L - C	428 - R - C	X Large	8½" +



### Carpal Tunnel Wrist Brace

- · Ideal for carpal tunnel syndrome, wrist sprains, strains or after cast removal
- · Maximum contouring at the palmer crease allows for full finger function
- · Removable and and adjustable palmer stay for maximum support
- Unique stockinette liner provides comfort and easy slip on application
- Velcro<sup>™</sup> closure for easy adjustment
- 8" length

Indications: Carpal tunnel syndrome, sprains, strains, tendonitis, wrist immobilization and post cast healing. | Suggested HCPC L3908

Order # Left	Order # Right	Size	Wrist Circumference
131 - L	121 - R	X Small	4½" - 5½"
133 - L	123 - R	Small	5½" - 6½"
135 - L	125 - R	Medium	6½" - 7½"
137 - L	127 - R	Large	7½" - 8½"
138 - L	128 - R	X Large	8½" +







#### **Universal Wrist Brace**

- Universal sizing
- · Adjustable dorsal stay
- · Constructed of lightweight breathable material

Indications: Carpal tunnel syndrome, sprains, strains, tendinitis, wrist wrist immobilization and post cast healing. | Suggested HCPC L3908

Order #	Size
340 - R	Right Universal
350 - I	Left Universal









### K.S. Lace-Up Wrist Brace

- · Deluxe ultra-suede perforated material with flannel lining for added patient comfort
- Two metal dorsal stays on the back of the brace helps prevent extension
- · Unique stockinette liner provides comfort and easy slip-on application

Indications: Carpal tunnel syndrome, sprains, strains, tendonitis, wrist immobilization and post cast healing. | Suggested HCPC L3908

Order # Left	Order # Right	Size	Wrist Circumference
351 - L	341 - R	X Small	4½" - 5½"
353 - L	343 - R	Small	5½" - 6½"
355 - L	345 - R	Medium	6½" - 7½"
357 - L	347 - R	Large	7½" - 8½"
358 - L	348 - R	X Large	8½" +







### Cock-Up Elastic Wrist Brace

- · Ultra lightweight style wrist brace
- · Removable and adjustable palmar stay for maximum support
- Velcro<sup>™</sup> closure for easy adjustment
- 7" length

Indications: Carpal tunnel syndrome, sprains, strains, tendonitis, wrist immobilization and post cast healing. | Suggested HCPC L3908

Order # Left	Order # Right	Size	Wrist Circumference
141 - L	111 - R	X Small	4½" - 5½"
143 - L	113 - R	Small	5½" - 6½"
145 - L	115 - R	Medium	6½" - 7½"
147 - L	117 - R	Large	7½" - 8½"
148 - L	118 - R	X Large	8½" +









### AO Easy Adjustable Wrist Orthosis®

- · Strong flexion synergy of the fingers and the wrist
- · Adjustable hinge at the wrist to help increase wrist extension
- · Optional finger separators help to prevent pressure areas
- · Fits either left or right wrist

Indications: Ideal for patients with tendinitis, arthritis, mild sprains and strains of the wrist. | Suggested HCPC L3916

Order # | Size

4100-E Universal



#### Wrist Hand Orthosis

- Designed for the treatment of severe strains and sprains to the wrist or early cast removal for fractures
- Cast-like rigidity provides ideal immobilization without the discomfort associated with traditional casts
- Foam tongue insert ensures comfort and additional padding for the styloid process
- · Complete with a double thickness cotton stockinette
- 8" length

Indications: Carpal tunnel syndrome, sprains, strains, tendonitis, wrist immobilization and post cast healing. | Suggested HCPC L3984

Order # Left	Order # Right	Size	Palm Circumferen	се
71063	71073	Small	2½" - 3"	回線線回
71065	71075	Medium	3" - 3½"	
71067	71077	Large	3½" - 4"	



#### Neoprene Wrist-Wrap Support

- · Neoprene coated with 4-way stretch nylon fabric
- · Provides compression and warmth to wrist area
- Velco<sup>™</sup> closure for easy adjustment
- · Fits either left or right wrist

Indications: Mild wrist sprains and strains.

Order # | Size | Universal







#### Elastic Slip-On Wrist Support

- · Durable knitted elastic provides support with flexibility
- · Easy slip-on design
- · Seamless construction will not irritate skin
- · Fits either left or right wrist

Indications: Mild sprains, strains, tendonitis, arthritis.

Order #	Size	Wrist Circum	ference	
1903	Small	4" - 6"		
1905	Medium	6" - 8"	and Report	具数线具
1907	Large	8" - 10"	1	
1908	X Large	10" - 12"	8 sq.,	回数操

#### **Elastic Wrist Guard Support**

- Unique 3" wide elastic material adheres to Velcro™ strap for exact fit and easy adjustment
- · Fits either left or right wrist

Indications: Mild sprains, strains, tendonitis, arthritis.

Order #	Size		
2900	Universal	D 447	

#### Thumb Spica Wrist Brace

- · Constructed of lightweight, breathable, padded material for maximum comfort and wear
- · Maximum contouring of the palmer crease allows for full finger function
- · Removable and adjustable palmar stay for maximum support
- · Removable rigid thumb stay is adjustable for a custom fit
- 8" length

Indications: Gamekeeper thumb, indicated for De Quervain's Syndrome and rheumatoid arthritis. | Suggested HCPC L3807 and L3809

Order # Left	Order # Right	Size	Wrist Circumference	
181-L	171-R	X Small	4½" - 5½"	
183-L	173-R	Small	5½" - 6½"	
185-L	175-R	Medium	6½" - 7½"	具数器具
187-L	177-R	Large	7½" - 8½"	
188-L	178-R	X Large	8½" +	





### Lycra Lined Wrist Brace with Thumb Spica

- · Soft, lycra lined interior for exceptional comfort
- · Lycra lining to wick away skin moisture
- · Removable palmer stay to allow correct fit

Indications: De Quervain's Syndrome, wrist and thumb sprains and strains, post cast healing. | Suggested HCPC L3807 and L3809

Order # Left	Order # Right	Size	Wrist Circumference
181 - L - C	171 - R - C	X Small	4½" - 5½"
183 - L - C	173 - R - C	Small	5½" - 6½"
185 - L - C	175 - R - C	Medium	6½" - 7½" <b>- 14333</b>
187 - L - C	177 - R - C	Large	7½" - 8½"
188 - L - C	178 - R - C	X Large	8½"+







### Universal Wrist Brace with Thumb Spica

- · Constructed of durable, lightweight material
- Universal sizing
- Adjustable dorsal stay
- · Removable and malleable aluminum palmer stay for correct fitting

Indications: Wrist immobilization, De Quervain's Syndrome, wrist and thumb sprains and strains. | Suggested HCPC L3807 and L3809

Order #

Size

170-R Right Universal

180-L

Left Universal





#### Premium Thumb Brace

- · Perforated suede thumb splint provides thumb immobilization and protection
- Constructed of lightweight, breathable, padded material for maximum comfort and wear
- · Perforated to allow moisture to escape the skin area
- · Malleable, rigid thumb stay is adjustable for a custom fit
- Velcro<sup>™</sup> closure for easy adjustment
- · Fits either left or right thumb

**Indications:** De Quervain's Syndrome, thumb sprains and strains. Suggested HCPC L3923

Order # | Size

21003 Universal





#### Universal Wrist / Forearm Brace

- Brace can be removed and put back on easily once properly fitted Velcro™ closure for easy removal
- 11" length

Indications: Tendonitis, arthritis, mild sprains and strains | Suggested HCPC L3908

Order #

SIZE

150 - R Right Universal

160 - L Left Universal



### Neoprene Wrist/Thumb Wrap Support

- · Neoprene coated with 4-way stretch nylon fabric
- · Provides compression and warmth to wrist/thumb area
- · Extended thumb provides abduction
- Velcro<sup>™</sup> closure for easy adjustment
- · Fits either left or right wrist/thumb

Indications: Mild wrist and thumb sprains and strains.

Order # 21002

012

Universal

















### Thumb Spica Orthosis

- Lightweight orthosis is designed to comfortably restrict motion of the basilar joint while allowing unobstructed motion of the fingers and wrists
- · May be modified using cast scissors and a heat gun for customized fit
- Complete with a double thickness cotton stockinette

Indications: Thumb injuries (ligament instability, sprain or muscle strain), gamekeeper's thumb, osteoarthritis, De Quervain's Syndrome, scaphoid fracture, thumb sprains. I Suggested HCPC L3807 and L3809

fracture, thumb	sprains.	Sugg	ested HCF	PC L3807 and L3809	
Order # Left	Order #	Right	Size	Palm Width	

1263	71253	Small	2½" - 3"	回線
1265	71255	Medium	3" - 31/2"	13/8
1267	71257	Large	3½" - 4½"	



#### Wrist Hand Thumb Orthosis

- Durable orthosis comfortably restricts motion of the wrist, hand and thumb
- Specially formulated thermoplaste provides cast-like rigidity without the discomfort associated with traditional casting methods

Indications: Indicated for use in post cast fracture removal, severe strains and sprains. | Suggested HCPC L3807 and L3809

Order # Left	Order # Right	Size	Palm Width
72063	72073	Small	2½" - 3"
72065	72075	Medium	3" - 3½"
72067	72077	Large	3½" - 4"



# 7 | 72077 | Large | 3½" - 4" 🔲 🔀

#### Thermoskin Arthritic Gloves

- · Provides warmth and equal compression to help relieve arthritis pain
- Velcro<sup>™</sup> closure allows added comfort and fit
- The outer layer is textured to provide additional grip

Indications: Pain and strains due to arthritis.

Order #	Size	Hand Circumference
82199	X Small	6" - 6¾"
83199	Small	7" - 7¾"
84199	Medium	8" - 9"
85199	Large	9¼" - 10½"
86199	X Large	10¾" - 11½"
87199	2X Large	11¾" +









#### Hinge R.O.M. Elbow Brace

- Telescopic length adjustments at forearm and biceps
- Plush padding for comfort fit
- Customizing arm and bicep cuffs for added comfort
- · Adjustable strapping system

Indications: Soft tissue strains / repairs, Post-op surgery. Suggested HCPC L3760

Order#

Size

2320-R Right Universal

2330-L Left Universal







#### Air Tennis Elbow System

- · Air tennis elbow system has three uniquely designed features
- · Air therapy pad helps reduce swelling
- · Foam tension pad to maximize compression and dissipate force
- · Fits either left or right forearm

Indications: Tennis elbow syndrome, tendinitis, lateral / medial epicondylitics.

Order # | Size

Universal







#### **Premium Tennis Elbow Support**

- · Constructed of 2" web strap
- Velcro<sup>™</sup> closure for easy application
- · Fits either left or right forearm

Indications: Tennis elbow syndrome, tendinitis, lateral/medial epicondylitics.

Order # |

Size

2300 Universal





### Band-IT™ Tennis Elbow Support

- · Provides intermittent pressure across the extension & supinator muscles
- Velcro<sup>™</sup> closure
- · Fits either left or right forearm

Indications: Tennis elbow syndrome, tendinitis, lateral/medial epicondylitics.

Order#

I Size

2302 Universal







#### Neoprene Tennis Elbow Sleeve with Strap

- Circumferential 2" Velcro<sup>™</sup> strap
- · Provides warmth and equal compression
- Velcro<sup>™</sup> closure
- · Fits either left or right forearm

**Indications:** Tennis elbow syndrome, tendonitis, lateral/medial epicondylitis.

Order #	Size	Forearm Circumference
2303	Small	9" - 10"
2305	Medium	10" - 11"
2307	Large	11" - 12"
2308	X Large	12" - 14

Measure the largest circumference of the forearm.



#### Elastic Slip-On Elbow Support

- · Easy slip-on
- · Seamless construction will not irritate skin
- · Provides intermediate compression to the elbow area
- · Durable knitted elastic provides support with flexibility

Indications: Elbow sprains and strains, tendonitis, bursitis.

Order #	Size	Forearm Circumference
2313	Small	8½" - 10"
2315	Medium	10" - 11½"
2317	Large	11½" - 13"
2318	X Large	13" - 14½"

Measure the largest circumference of the forearm.



#### Elbow and Heel Protector

- · Most popular, effective elbow and heel protector
- · Protector helps reduce pressure on sensitive or injured area
- · Lightweight mesh tubular sleeve
- No seams in critical areas
- · Cool padded inner cushion prevents "bottoming out"
- · Contours extremely well to elbow or heel for added protection
- · Packaged one pair per box

Indications: Sore, tender elbow injuries that require padded protection.









### Aero Walker® (High Top and Low Top)

- Inflatable bladder incorporated into liner for increased immobilization of the lower leg, ankle, and foot
- Reinforced adjustable toe guard for added protection
- Cushioned sole for greater energy absorption
- Mold injected rocker bottom promotes a durable natural gait pattern
- Fits either left or right foot

Indications: Used for soft tissue injuries, grade 2 and 3 sprains, stable fractures and post-operative stabilization. | Suggested HCPC L4360 and L4361

Order # High Top	Order # Low Top	Size	Male Shoe Size	Female Shoe S	Size
320 - AG	320 - LAG	X Small	Up to 4	Up to 5½	
330 - AG	330 - LAG	Small	4½ - 7	6 - 8	
360 - AG	360 - LAG	Medium	7½ - 10	8½ - 11½	15 <del>88</del> 8
390 - AG	390 - LAG	Large	10½ - 12	11½ - 13½ 🔻	100
398 - AG	398 - LAG	X Large	12½ +	13½ +	





- · Lightweight for patient comfort
- · Cushioned heel for energy absorption
- · Anatomically correct design offers superior patient comfort
- Rocker bottom design for natural gait pattern
- · Molded plastic uprights
- Mold injected rocker bottom promotes a durable natural gait pattern
- · Fits either left or right foot

Indications: Used for soft tissue injuries, grade 2 and 3 sprains, stable fractures and post-operative stabilization. | Suggested HCPC L4386 and L4387

			00	
Order # High Top	Order # Low Top	Size	Male Shoe Size	Female Shoe Siz
320 - Z	320 - LZ	X Small	Up to 4	Up to 51/2
330 - Z	330 - LZ	Small	4½ - 7	6 - 8
360 - Z	360 - LZ	Medium	7½ - 10	8½ - 11½
390 - Z	390 - LZ	Large	10½ - 12	11½ - 13½
398 - Z	398 - LZ	X Large	12½ +	13½ +





- Inflatable bladder incorporated into liner for increased immobilization of the lower leg, ankle and foot
- · Lightweight, low profile design for greater patient comfort
- · Cushioned heel for greater energy absorption
- · Mold injected rocker bottom promotes a durable natural gait pattern
- Fits either left or right foot and Pink Liner available for Low Top models\*

Indications: Used for soft tissue injuries, grade 2 and 3 sprains, stable fractures and post-operative stabilization. Helps to reduce pain, swelling and edema of foot and ankle. | Suggested HCPC L4360 and L4361

Order # High Top	Order # Low Top	Size	Male Shoe Size	Female Shoe Size
320 - AZ	320 - LAZ	X Small	Up to 4	Up to 5½
330 - AZ	330 - LAZ	Small	4½ - 7	6 - 8
360 - AZ	360 - LAZ	Medium	7½ - 10	8½ - 11½ <b>- 222</b>
390 - AZ	390 - LAZ	Large	10½ - 12	11½ - 13½
398 - AZ	398 - LAZ	X Large	12½ +	13½ +









#### **Pediatric Walker**

- · Excellent for infants
- · Extremely lightweight
- Soft, plush lining
- Fits left or right foot

Indications: Used for soft tissue injuries, grade 2 and 3 sprains, stable fractures and post-operative stabilization. | Suggested HCPC L4386 and L4387

Order #	Size	Years of Age	Sole Length	
330 - P	Small / Medium	1 - 2	5½"	回線深间
390 - P	Large	2½ - 6	61/4"	
398 - P	X Large	6½ - 9½	71/4"	



### Low Profile Air Walker (High Top and Low Top)

- · Metal uprights
- Inflatable bladder incorporated into liner for increased immobilization of the lower leg, ankle and foot
- · Lightweight, low profile design for greater patient comfort
- Fits either left of right foot

Indications: Used for soft tissue injuries, grade 2 and 3 sprains, stable fractures and post-operative stabilization. Helps to reduce pain, swelling and edema of foot and ankle. | Suggested HCPC L4360 and L4361

Order # High Top	Order # Low Top	Size	Male Shoe Size	Female Shoe Size
330 - AE	330 - LAE	Small	4 - 6	5 - 7
360 - AE	360 - LAE	Medium	7 - 10	8 - 11
390 - AE	390 - LAE	Large	11 - 13	12 - 14



### Low Profile Walker (High Top and Low Top)

- · Metal uprights
- · Lightweight for patient comfort
- · Rocker design promotes a smooth energy conserving gait
- Cushioned heel for energy absorption
- · Anatomically correct design offers superior patient comfort
- Fits either left or right foot

**Indications:** Used for soft tissue injuries, grade 2 and 3 sprains, stable fractures and post-operative stabilization.

Suggested HCPC L4386 and L4387

Order # Hi Top	Order # Low Top	Size	Male Shoe Size	Female Shoe Size
330 - E	330 - LE	Small	4 - 6	5 - 7
360 - E	360 - LE	Medium	7 - 10	8 - 11
390 - E	390 - LE	Large	11 - 13	12 - 14



### Stabilizer Range of Motion Walker

- · Patented walker design offers a superior rocker bottom in a low profile walker
- · Continuous arching rocker design promotes a smooth energy conserving gait
- Adjustable ankle motion ranges from fixed to 45°
- · Fits either left or right foot

Indications: Used for soft tissue injuries, grade 2 and 3 sprains, stable fractures and post-operative stabilization. | Suggested HCPC L4386 and L4387

Order #	Size	Male Shoe Size	Female Shoe Size
313	Small	4 - 7	5 - 81/2
315	Medium	7½ - 10	9 - 11
317	Large	10½ - 12	11½ - 13
318	X Large	12½ +	13½ +

#### P.F.S. Night Splint

- Allows the patient to apply consistent stretch to the plantar fascia by holding the foot in gentle dorsiflexion
- · Foam wedge is included for additional stretch on the plantar fascia
- · Lightweight design for patient comfort while sleeping
- · Fits left or right foot

Indications: Designed to treat plantar fasciitis, Achilles tendonitis and heel pain. | Suggested HCPC L4396 and L4397

Order #	Size	Male Shoe Size	Female Shoe Size
4833	Small	4 - 6	5-7
4855	Medium	6 - 10	7 - 11
4877	Large	10 - 13	11 - 14

#### **Dorsal Night Splint**

- · Flexible, breathable material provides a comfortable fit
- · Very well padded
- · Fits right or left foot

**Indications:** Designed to treat plantar fasciitis, Achilles tendonitis and heel pain. | Suggested HCPC L4396 and L4397

			_
Order #	Size		Female Shoe Size
	Small / Medium		Up to 10
4978	Large / X Large	91/2+	101/4+





Updated

Design

#### **Podus Boot**

- Helps promote healing of foot and ankle pressure sores
- · Helps prevent foot drop and helps eliminate heel pressure
- Helps limit hip and leg rotation
- · Fits either left or right foot

Size

Non-ambulatory

Order#

**Indications:** For the treatment and prevention of ankle/foot contractures and internal/external rotation of the lower extremity. | Suggested HCPC L4396 and L4397

3300	Universal (8" - 13" approx till size 16 in Male Shoes Size)
	*Calf Circumference: Un to 16"













## AO Carbon Fiber AFO

- At toe off, the limb receives energy and propulsive help
- · During swing phase, the foot is supported for ground clearance
- · Carbon fiber design offers a high strength-to-weight ratio as well as dynamic mobility
- · At first contact, plantar flexion is gradual
- · Ground response forces imparted to the foot correspond to the stiffness gradient at the heel and toe

Indications: Drop foot secondary to: - Cerebrovascular accident (CVA) - other neurological conditions - multiple sclerosis - mild knee instability - Charcot Marie Tooth Disease (HMSN) overuse injuries to the lower extremities.

#### Posterior/Lateral, L1951:

1	Order # Left	Order # Right	Size	Male Shoe Size	Female Shoe Size
		7511	X-Small	Up to 2 ½ 3 - 6	Up to 4 ½ 5 - 8
	7523	7513	Small	3 - 6	
	7525	7515	Iviedium	6-9	8 - 11
		7517	Large	9 - 11	11 - 13

#### Posterior/Medial, L1951:

ı	Order # Left	Order # Right	Size	Male Shoe Size	Female Shoe Size
	7623		Small	3 - 6	5 - 8 8 - 11
	7625	7615	Medium		
			Large	9 - 11	11 - 13

#### Anterior/Lateral, L1932:

Order # Left	Order # Right	Size	Male Shoe Size	Female Shoe Size
	7711	X-Small	Up to 2 ½ 3 - 6 6 - 9	Up to 4 ½ 5 - 8
7723	7713	Small	3 - 6	5 - 8
7725	7715 7717	Medium	6 - 9	8 - 11 11 - 13
7727	7717	Large	9 - 11	11 - 13







# Ankle Foot Orthosis (AFO)

- 13" height
- · Provides full range of plantar flexion and dorsiflexion
- Trimmed low in the malleolus area
- · Designed to prevent foot drop
- · Designed to fit into the shoe

Indications: Designed to aid in foot drop that is secondary to CVA or nerve palsy. | Suggested HCPC L1930

Order #   Left	Order #   Right	Size	Male Shoe Size	Female Shoe	Size
7023	7013	Small	3 - 6	5 - 8	
7025	7015	Medium	6 - 9	8 - 11	回胡絮
7027	7017	Large	9 - 11	11 - 13	
7028	7018	X Large	11+	N/A	







- Prevents abnormal ankle inversion, eversion
- · Offers more stability than standard ASO ankle braces
- · Low profile
- · Light weight

Indications: Ankle injury rehabilitation, Sub-acute ankle sprains, chronic ankle instabilities or prophylactic use.

Suggested HCPC L1971 (with Attached Calf Extension)

Suggested HCPC L1906 (without Attached Calf Extension)

Order # Left	Order # Right	Size	Male Shoe Size	Female Shoe Size
833 - L	823 - R	Small	Up to 8	Up to 9
835 - L	825 - R	Medium	8½ - 12	9½ - 13
837 - L	827 - R	Large	12+	13+







## Air-Gel Ankle Brace

- · The Air-Gel Brace with its pneumatic and gel chambers are state-of-the art
- Removable Air-Gel bladders can provide for "cold" therapy
- · Helps eliminate swelling and edema
- For acute injuries when ankle is still tender or swollen
- · Use when either walking or jogging
- Fits either left or right ankle

**Indications:** Sub-acute ankle sprains, chronic ankle instabilities or prophylactic use. | Suggested HCPC L4350

Order #	Size	Brace Height	回数数数
442	Regular	10½"	
442 - C/T	Child / Trainer	9"	



## Air-Lite Ankle Brace

- · Medical grade air bladder
- Air-Lite Brace may be inflated to patients comfort
- · Contoured shells conform closely to ankle
- Athletic performance
- · Fits either left or right ankle

**Indications:** Sub-acute ankle sprains, chronic ankle instabilities or prophylactic use. | Suggested HCPC L4350

Order #	Size	Brace Height		
443-P	Universal	10%"	BAST	



## Cool Foam Ankle Brace™

- Ultra lightweight functional
- Memory foam interface padding provides maximum comfort
- Designed to help support and protect unstable or injured ankles
- · Allows normal dorsal and plantar flexion while controlling inversion and eversion
- Fits either left or right ankle

**Indications:** Sub-acute ankle sprains, chronic ankle instabilities or prophylactic use. | Suggested HCPC L4350

Order #	Size	Brace Height	■:
441	Regular	10½"	- 34
441 - C	Child	81/2"	





## Quick Lace Ankle Brace

- Lightweight and low profile design that is highly breathable with quick drying fabric
- Top strap design for increased stabilization around shin
- · Speed laced design provides the option of a quick lace-up
- Figure 8 finger pull loops allow the user to easily adjust the strapping during donning and for a quick adjustment during use
- · Fits either left or right ankle

Indications: Chronic ankle instability, sub-acute ankle sprains, ankle injury rehabilitation. | Suggested HCPC L1902

Order #	Size	Male Shoe Size	Female Shoe Size	
361	X Small	6 - 7	7 - 8	
363	Small	7 - 9	8 - 10	
365	Medium	9 - 11	10 - 12	₽
367	Large	11 - 13	12 - 14	,
368	X Large	13 - 14	14 - 15	Ö



## **AO Stabilizer Ankle Brace**

- Advanced ankle design provides the compression and comfort of a soft ankle brace but with 20% - 35% more protection than traditional stirrups
- Lightweight and low profile design with highly breathable and quick drying fabric
- · Speed laced design provides the option of a quick lace-up
- Figure 8 finger pull loops allow the user to easily adjust the strapping during donning and for a quick adjustment during use
- Fits either left or right ankle

**Indications:** Chronic ankle instability, sub-acute ankle sprains, ankle injury rehabilitation. | Suggested HCPC L1902

Order #	Size	Male Shoe Size	Female Shoe	e Size	
373	Small	6 - 7	7 - 8	_	
375	Medium	7 - 10	8 - 11	de Broce	
377	Large	10 - 12	11 - 13		
378	X Large	12 - 14	13 - 15	9 332	186762



## Advanced Ankle Helper Hinge Brace

- Modified orthotic footplate provides critical foot and arch support, resisting inversion and eversion
- · Neutral arch position promotes healing of injured ligaments
- Three strap adjustment system provides custom fit

Indications: Ankle injury rehabilitation, Sub-acute ankle sprains, chronic ankle instabilities or prophylactic u e. Mild drop foot. | Suggested HCPC L1906

Order # Left	Order # Right	Size	Male Shoe Size	Female Shoe Size
19-0020L	19-0021R	X Small	N/A	4 - 6
19-0030L	19-0031R	Small	4 - 7	7-9
19-0040L	19-0041R	Medium	8 - 10	10 - 12
19-0050L	19-0051R	Large	11 - 13	13 - 15





## Silicone Elastic Ankle Support

- Viscoelastic silicone insert helps provide relief from pain, swelling and inflammation
- · Three dimensional knotted elastic
- 7" sleeve

**Indications:** Sub-acute ankle sprains, chronic ankle instabilities or prophylactic use. | Suggested HCPC L1901

Order #	Size	Ankle Circumference	
453	Small	71/8" - 83/8"	
455	Medium	83/8" - 91/4"	
457	Large	91/4" - 10"	
458	X Large	10" - 10¾"	



# Elastic Slip-On Ankle Support

- Durable lightweight knitted elastic provides support with flexibility
- · Unique seamless construction will not irritate skin
- · Open heel and toe provides easy slip on
- · Fits either left or right ankle

Indications: Sub-acute ankle sprains, chronic ankle instabilities or prophylactic use. | Suggested HCPC L1901

Order #	Size	Ankle Circum	nference	
2703	Small	7" - 8"		
2705	Medium	8" - 9"	Seeks and the	■総総■
2707	Large	9" - 10"		
2708	X Large	10" - 11"	6 sag-	



## Thermoskin Elastic Ankle

- · Anatomically designed with no seam
- · Elastic support will comfortably contour to the shape of your ankle
- Provides equal compression and support for weak and injured ankles, while minimizing bunching

Indications: Sub-acute ankle sprains, chronic ankle instabilities or prophylactic use. | Suggested HCPC L1901

Order #	Size	Ankle Circumference	
83604 - R	Small	61/4" - 81/4"	
84604 - R	Medium	8¼" - 10¼"	回数线回
85604 - R	Large	10¼" - 12¼"	
86604 - R	X Large	12¼" - 14¼"	回路線



## **Full Insole Silicone Foot Orthosis**

- Covers the full sole and helps prevent fatigue due to its massaging effect which improves blood circulation
- · Designed to relieve pressure on muscles, tendons and joints
- Colored striking areas are softer than surrounding transparent material giving flexible pressure relief
- · Constructed of 100% silicone

Indications: Heel spurs, metatarsal pain, calluses and fat pad atrophy.

Order #	Size	Male Shoe Size	Female Shoe Size
053	Small	5 - 6	6 - 7
055	Medium	6 - 7	7 - 8
057	Large	7 - 9	8 - 10
058	X Large	9 - 11	10 - 12
059	2X Large	11 - 13	



## Three Quarter Silicone Foot Orthosis

- Covers three quarters of the sole & helps prevent fatigue due its massaging effect which improves blood circulation
- · Designed to relieve pressure on muscles, tendons and joints
- Colored striking areas are softer than surrounding transparent material, giving flexible pressure relief
- Constructed of 100% silicone

**Indications:** Heel spurs, metatarsal pain, calluses and fat pad atrophy.

Order #	Size	Male Shoe Size	Female Shoe Size
063	Small	5 - 8	6 - 9
065	Medium	8 - 10	9 - 11
067	Large	10 - 13	11 - 14



# **Central Heel Support**

- · Helps relieve lower extremity fatigue through shock absorption of the heel area
- · Designed to relieve pressure on muscles, tendons and joints
- Colored striking areas are softer than surrounding transparent material giving flexible pressure relief
- · Constructed of 100% silicone

**Indications:** Plantar fasciitis, Achilles tendonitis, calcaneal and heel spurs. | Suggested HCPC L3170

Order #	Size	Male Shoe Size	Female Shoe Size
083-E	Small Medium	5 - 8	6 - 9
085-E	Medium	8 - 10	9 - 11
087-E	Large	10 - 13	11 - 14





## Diabetic Insert

 Dual density diabetic insert is heat moldable at temperature of 230°F or higher

Indications: Full length foot orthotic for patients at risk to skin breakdown and ulceration due to diabetes or decreased sensation in extremities. | Suggested HCPC A5512

Order #	Size	Male Shoe Size	Female Shoe Size
471	XSmall	N/A	5 - 7
473	Small	6 - 7	8 - 9
475	Medium	8 - 9	10 - 11
477	Large	10 - 11	12 - 13
478	X Large	12 - 13	N/A
479	2X Large	14 - 15	N/A



# Tuli Heel Cup

- New double-ribbed heel cups provide 50% more protection
- · Helps relieve arch, ankle, leg and back pain
- Packaged as one pair per box

Indications: Relief from Plantar Fasciitis heal pain and shin splints.

ı	Order #	Size
	100 - H	Regular - Under 175 lbs.
	101 - H	Large - Over 175 lbs.



## Elbow and Heel Protector

- · Most popular, effective elbow and heel protector
- · Protector helps reduce pressure on sensitive or injured area
- · Lightweight mesh tubular sleeve
- · No seams in critical areas
- Cool padded inner cushion prevents "bottoming out"
- · Contours extremely well to elbow or heel for added protection
- Packaged as one pair per box

**Indications:** To prevent bruising, sores and ulcerations. Not suitable to cover open wounds.

Order #	Size
2310	Universal





# Classic Post-Op Shoe

- · Flared, semi-rigid out sole provides medial and lateral stability
- Breathable nylon mesh upper
- Simple, one buckle closure on each strap

#### Suggested HCPC L3260

33			
Order #	Size   Unisex	Male Shoe Size	Female Shoe Size
CL1	Pediatric	10 - 12	10 - 12
CL2	Women's Small	N/A	4 - 6
CL3	Women's Medium	N/A	6½ - 8
CL4	Women's Large	N/A	8½ - 10
CL5	Men's Small	6 - 8	N/A
CL6	Men's Medium	8½ - 10	10½ -12
CL7	Men's Large	10½ - 12	12½ - 14
CL8	Men's X Large	12½ - 14	N/A 🌄
CL9	Men's 2X Large	14 - 16	N/A



# Original Med-Surg<sup>™</sup> Shoe

Order # Male	Order # Female	Size	Male Shoe Size	Female Shoe Size
MSM1N	MSW1N	Small	6 - 8	4 - 6
MSM2N	MSW2N	Medium	8½ - 10	6½ - 8
MSM3N	MSW3N	Large	10½ - 12	8½ - 10
MSM4N	N/A	X Large	12½ - 14	N/A





# Open Toe Cast Boot

Order #	Size   Unisex	Male Shoe Size	Female Shoe Size
CB00	X Small	N/A	N/A
CB01	Small	N/A	5½ - 7½
CB02	Medium	6 - 81/2	8 - 10
CB03	Large	9 - 11	N/A
CB04	X Large	12 - 14	N/A

Suggested HCPC L3260



# Ortho Wedge™ Shoe

9				
Order #	Size   Unisex	Male Shoe Size	Female Shoe Size	
OW1B	Small	6 - 8	7½ - 10	
OW2B	Medium	8½ - 10	10½ - 13	
OW3B	Large	10½ - 12	N/A	
OW4B	X Large	12 ½ - 14	N/A	

Suggested HCPC A9283



# Heel Wedge™ Shoe

Order #	Size   Unisex	Male Shoe Size	Female Shoe Size
HW1B	Small	6 - 8	7½ - 10
HW2B	Medium	8½ - 10	10½ - 13
HW3B	Large	10½ - 12	N/A
HW4B	X Large	12½ - 14	N/A

Suggested HCPC A9283





# Deluxe Universal Length Arm Sling

- Envelope style
- Can be folded to proper length to accommodate most patients
- Soft foam neck pad for greater patient comfort
- · Slide buckle on strap for easy adjustment
- Velcro<sup>™</sup> closure

Indications: Arm injuries or instabilities. Post surgery.

Order # | Size

2200 Universal







# Velpeau Immobilizer

- · Breathable cool cotton material
- Webbing swathe strap
- Soft foam neck pad for greater patient comfort
- · Slide buckle on strap for easy adjustment
- Velcro<sup>™</sup> closure

**Indications:** Shoulder injuries or instabilities. Post Surgery. Suggested HCPC L3660

Order #	Size	Strap Length
2203	Small	36"
2205	Medium	39"
2207	Large	43"
2208	X Large	46"



# **Premium Arm Sling**

- · Our most popular arm sling
- Velcro<sup>™</sup> closure
- · Easy to adjust

Indications: Arm injuries or instabilities post surgery.

Order #	Size		
2211	X Small		
2213	Small		
2215	Medium	Am Sing	国统统
2217	Large		
2218	X Large	Ø 800°	
	2211 2213 2215 2217	2213 Small 2215 Medium 2217 Large	2211 X Small 2213 Small 2215 Medium 2217 Large





## Shoulder Pillow with Ball

- · Special design fits right or left arm
- High density foam pillow with attached stress ball helps promote hand circulation without moving arm or shoulder

Indications: Glenohumeral dislocations/subluxation, soft tissue repair, posterior shoulder repair. | Suggested HCPC L3670

Order #	Size	Length from wrist to e	lbow
2903	Small	9" - 12"	
2905	Medi	12" - 15"	
2907	Large	15" - 18"	
2908	X Large	18" - 21"	



# Deluxe Sling and Swathe

- · Cool, lightweight, pile/foam construction provides comfort
- Velcro<sup>™</sup> closure for easy adjustment
- Fits either right or left arm

Indications: Shoulder dislocations, injuries or instabilities.

Suggested HCPC L3670

Order #	Size	Swathe Length	
2913	Small	34"	
2915	Medium	40"	回数300回
2917	Large	48"	
2918	X Large	54"	



# **Deluxe Shoulder Abduction System**

- · Lightweight, universally-sized
- Removable wedge allows shoulder positioning at 0, 15, 30, 60 and 90 degrees
- Strapping systems allows abductor to be worn on the left or right arm with minimal adjustment

Indications: Post-operative rotator cuff, elbow and shoulder injuries. Rotator cuff tear and tendinitis at the shoulder. Adhesive capsulitis at the shoulder (frozen shoulder). Cumulative trauma/repetitive strain injury. Elbow fractures including radial head. Stroke and CVA recovery and residual hemiparesis. | Suggested HCPC L3960

Order # | Size 5100-E | Universal





## **Elastic Shoulder Immobilizer**

- · Fits left or right arm
- Elastic waist band
- · Easy adjustment, lightweight and comfortable

**Indications:** Shoulder injuries, instabilities or dislocations. Suggested HCPC L3650

Order #	Size	Rib Circumference	
2811	X Small	23" - 27"	
2813	Small	27" - 31"	
2815	Medium	31" - 35"	٥
2817	Large	35" - 39"	١
2818	X Large	39" - 43"	Į



## Universal Shoulder Immobilizer

- · Immobilizes arm and shoulder
- · Contoured chest strap fits both male and female
- Shoulder strap and forearm cuff may be placed in desired position
- $^{ullet}$  Velco  $^{\text{TM}}$  closure for easy application
- Fits either left or right arm

Indications: Shoulder dislocations, injuries or instabilities.
Suggested HCPC L3650

Order #	Size	
2800	Universal	



# **Premium Clavicle Support**

- Our most popular Velcro<sup>™</sup> closure clavicle support
- Comfortable wide padding for proper immobilization
- · Felt back pad helps prevent irritation
- Velcro<sup>™</sup> closure for easy application

**Indications:** Clavicle fracture, postural control. Suggested HCPC L3650

Order #	Size	Description	Chest Circumference
2603	Small	Child / Youth	25" - 30"
2605	Medium	Youth / Small Adult	30" - 35" 🚟 🖳
2607	Large	Adult	36" - 40"
2608	X Large	Large Adult	40" - 45"





# Therapeutic Aides



# Contoured **Back Rest**

- · Helps maintain correct sitting posture
- · Elastic strapping system fits almost any seat

Indications: Provides firm support to the lower region of the spine.

Order # | Size 14" x 15" 554



## Double Lobe **Cervical Pillow**

- · Supports the neck in a
- natural position
- Removable, washable cover

Indications: Provides support and pain relief to cervical area.

Order # | Size 16" x 22" 556



#### **Cervical Pillow**

- · Soft center cradles the
- · neck to help correct poor
- sleeping posture
- · Removable, washable

Indications: Provides support and pain relief to cervical area

Order #	Size
553	8" x 16"



# Aluminum Adjustable **Push-Button Crutches**

- Aluminum adjustable crutch with push-button height adjustment
- · Aluminum security nut and all accessories included
- 8 Pieces per case (16 pieces sold in case qty. only)

Order #	Size	Crutch Height	Patient Height
612	Youth	36" - 44"	4'6" - 5'2"
614	Adult	44" - 52"	5'2" - 5'10"
618	Tall / Adult	52" - 60"	5'10" - 6'6"

# Aluminum Adjustable **Folding Cane**

- · Carved wood handle
- Anodized black finish
- 300 lbs. weight limit

Order #	Size
650	33½" x 36½'



## Aluminum Adjustable Cane

- Attractive finish resists
- scratches
- 1" diameter anodized,
- extruded aluminum tubing

T.E.N.S Lead Wires

300 lbs. weight limit

Order #		Quantity
651	Adult (adjusts from 29" - 38")	6 / Case

## T.E.N.S. Unit

- · Dual Channel Unit
- · Adjustable pulse duration and frequency
- · Continuous burst and modulation setting
- · Includes carrying case, lead wires, electrodes and battery





#### Order # 3000-W-NS



## Twin Stim® TENS and EMS Combo

- · Five TENS Modes: Burst / Normal /
- Modulation / SD1 / SD2
- · Two EMS Modes: Synchronous,
- Alternate
- · System includes devices, lead wires, four self-adhesive reusable electrodes
- · 9-Volt battery, hard plastic carrying case and instruction manual
- Patient compliance meter can store up to 60 sets of operation manual



Order# 3001

# **Premium Self-Adhesive Electrodes**

- Comfortable, selfadhering electrodes
- Four per pack

Order #	Size
3014	2" X 2"







# **Dry Pro Cast Cover**

- · Most popular cast cover
- · Patented watertight vacuum seal
- · Made of high quality surgical latex
- · Non-skid-grid helps prevent slipping

Order # Arms	Size	Length	Circumference
HA - 13	Half Arm Small	17½"	7¾" - 10"
HA - 15	Half Arm Large	20"	10+"
FA - 12	X Small	16"	6" - 6¾"
FA - 14	Small	23"	6¾" - 8¾"
FA - 16	Medium	28"	8¾" - 10"
FA - 18	Large	31½"	10"+

Order # Legs	Size	Length	Circumference
HL - 13	Half Leg Small	21"	10" - 13"
HL - 15	Half Leg Large	23½"	13+"
FL - 12	X Small	19"	7¾" - 11"
FL - 14	Small	29"	14" - 16½"
FL - 16	Medium	33"	16½" - 21"
FL - 18	Large	37"	21"+



# Cast Protector & Bandage Protector

- Unique patented ring eliminates need for strapping or Velcro™ enclosure
- Latex free diaphragm stretches easily over cast to form watertight seal
- Durable polyvinyl protector prevents water penetration

Adult Sizes	Order #	Pediatric Sizes	Order #
Foot / Ankle (11")	2230	Medium Arm (18")	2105
Hand (11")	2229	Large Arm (22")	2107
Short Arm (22")	2223	Medium Leg (18")	2115
Long Arm (39")	2225	Large Leg (31")	2117
Short Leg (23")	2227		
Long Leg (42")	2228		
Extra Wide Opening	Order.#	Extra Wide Opening	Order.#
Short Arm (22")	6012	Short Leg (23")	2106

# Casting Tape

- Provides the benefits of a fiberglass casting tape with the handling ease of plaster to provide casting personnel with the casting characteristics they need
- 10 rolls per box

Order #	Description	
542	2" x 4" Yds	
543	3" x 4" Yds	
544	4" x 4" Yds	
545	5" x 4" Yds	







# Plaster Bandages

- Fast Set Time: 5-8 Minutes. (Blue Color)
- 12 Rolls per Box / 6 Boxes per case

Order # Description Quantity
41845 4" x 5" Yds 6 Doz. / Case



# **Cotton Stockinette**

- · Used under casts for skin protection
- Excellent memory retention
- Made from the softest, highest quality yarn available

Order #	Description
102	2" x 25" Yds
103	3" x 25" Yds
104	4" x 25" Yds
106	6" x 25" Yds



# **Elastic Bandages**

- · High-quality elastic bandage
- · Weave provides consistent compression
- Available with clips or Velcro closure

Order # w/Clips	Order # w/Velcro	Length	
3320-1	5925-1	2" X 5½" Yds	
3330-1	5935-1	3" X 5½" Yds	-
3340-1	5945-1	4" X 5½" Yds	900
3360-1	5965-1	6" X 51/2" Yds	Ø su:



<sup>\*10</sup> rolls per box w/o Retail Packaging (no -1 in the Order #)

<sup>\*1</sup> roll per box w/ Retail Packaging (with -1 in the Order #)

# **Aluminum Finger Splint Combination Kit**

- Kit contains nine frog splints, nine finger protector splints, nine baseball splints, nine curved finger splints and twelve finger cot splints
- 48 piece assortment kit

Order #	Size	
250	Combination	Ki



# **Curved Finger Splint**

- · Used to immobilize the joint following dislocation
- Used in the fractures of the phalanges where extension is required
- Each size packaged 12 per bag

Order #	Size	Length
213	Small	1½"
215	Medium	3"
217	Large	6"
218	X Large	9"



## **Finger Protector Splint**

- Used for injuries to the finger tips
- · Splints can be bent to any degree
- Each size packaged 12 per bag

Order #	Size	Length
233	Small	1½"
235	Medium	31/4"
237	Large	3½"





# Frog Splint

- · Maintains phalangeal joints in position
- · Malleable arms fold to hold finger in desired alignment
- Each size packaged 12 per bag

Order #	Size	Length
203	Small	21/8"
205	Medium	2¾"
207	Large	27/8"



# Finger Cot Splint

- · Used for protection to the finger tips
- Provides total finger immobilization
- · Each size packaged 12 per bag

Order #	Size	Length
221	X Small	1½"
223	Small	2"
225	Medium	3"
227	Large	3¾"



# **Baseball Finger Splint**

- · Stabilizes fractures of the distal phalanx
- Can be easily molded to give the desired amount of extension or hyper-extension
- Each size packaged 12 per bag

Order #	Size	Length
243	Small	3"
245	Medium	4"
247	Large	5"




## Foam Finger Strip

- Lightweight splint can be cut to size and conformed to the patient's needs
- Each size packaged 12 per box
- · Immobilize finger to correct position

Order #	Length
251	½" x 18"
252	¾" x 18"







## Men's Support Socks

- Ribbed Style
- Attractive Dress Style
- Comfortable Balloon Toe and Heel Pocket
- Thin, Lightweight Nylon/Spandex
- · Non-Binding, Comfort Top Band

#### 15-20 mm Hg Compression | 20-30 mm Hg Compression

Specify colors: BL-Black, BR-Brown, N-Navy, T-Tan Example: 9303-N

#### 30-40 mm Hg Compression

Specify colors: BL-Black, N-Navy Example: 9503-BL

Order #	Small	Medium	Large	X Large
15-20	9303	9305	9307	9308
20-30	9403	9405	9407	9408
30-40	9503	9505	9507	9508
Shoe Size	7 or less	7½ - 10	10½ - 12	12½ +
Ankle	7" - 9"	7½" - 9"	10" - 14"	11" - 15"
Calf	10½" - 15"	11" - 16½"	12" - 17"	13" - 19"
Length*	11" - 14"	12" - 15"	13" - 16"	14" - 18"



# Ladies' Support Socks

- Ribbed-Style Ladies Trouser Dress Socks
- Comfortable Balloon Toe and Heel Pocket
- Thin, Lightweight Nylon/Spandex
- · Non-Binding, Comfort Top Band

#### 15-20 mm Hg Compression

Specify colors: BL-Black, T-Tan, W-White Example: 9313-T

	Small	Medium	Large	X Large
Order #	9313	9315	9317	9318
Shoe Size	4 - 5	5½ - 7½	8 - 10½	11+
Ankle	61/4" - 71/4"	7½" - 9"	91/4" - 103/4"	11" - 11½"
Calf	11" - 13½"	13" - 15½"	15" - 17½"	17¾" - 18½
Lanath*	12" - 18"	13" - 10"	14" - 20"	15" - 22"



## Ladies' Knee Highs

- Sheer Therapy
- · Light, Silky Look and Feel
- · Reinforced Toe and Foot
- Thin, Lightweight Nylon / Spandex
- Non-Binding, Comfort Top Band

#### 15-20 mm Hg Compression

Specify colors: BL-Black, F-Fawn, NU-Nude

Example: 9323-BL

	Small	Medium	Large	X Large	2X Large
Order #	9323	9325	9327	9328	9329
Shoe Size	4 - 5	5½ - 7½	8 - 10½	11+	11+
Ankle	6¼" - 7¼"	7½" - 9"	9¼" - 10¾"	11" - 11½"	11¾ - 12¾
Calf	11" - 13½"	13" - 15½"	15" - 17½"	17¾" - 18½"	18 - 19½"
Length*	12" - 18"	13" - 19"	14" - 20"	15" - 22"	15 - 22"

<sup>\*</sup> Length is from heel to bend in knee.







<sup>\*</sup> Length is from heel to bend in knee.

<sup>\*</sup> Length is from heel to bend in knee.

# Ladies' Thigh Highs

- · Sheer Therapy
- · Lace Top with Silicone Beads
- · Reinforced Toe and Foot

#### 15-20 mm Hg Compression

Specify colors: BL-Black, F-Fawn, NU-Nude Example: 9343-NU

	Small	Medium	Large	X Large	2
Order #	9333	9335	9337	9338	ξ
Shoe Size	4 - 5	5½ - 7½	8 - 10½	11+	1
Ankle	6¼" - 7¼"	7½" - 9"	9¼" - 10¾"	11" - 11½"	1
Calf	11" - 13½"	13" - 15½"	15" - 17½"	17¾" - 18½"	1
Thigh	14" - 15½"	16" - 17½"	17½" - 19"	19" - 21"	2
Length	23" - 25"	25" - 28"	28" - 31"	31" - 34"	3





## Ladies' Pantyhose

- · Sheer Therapy
- Control Top
- · Reinforced Toe

#### 15-20 mm Hg Compression

Specify colors: BL-Black, F-Fawn

Specify	COIO15. DL-	Diack, F-F	avvii,
NU-Nud	e Example	: 9333-F	

Order #						100
Pantyhose	9343	9345	9347	9348	9349	-l



# **Anti-Embolism Stockings**

- · Available in Knee or Thigh High Lengths
- Thrombo-Embolic Deterrent
- · Use for Pre and Post-Operative
- · Available in Open or Closed Toe
- · Non-Binding, Comfort Top Band

#### 18 mm Hg Compression

Specify colors: W-White, BE-Beige, BL-Black Example: 9353-BE











# **Knee Highs**

- For Men and Women
- Available in Open or Closed Toe
- Comfortable Top Band
- Traditional, Thicker, Non-Sheer Fabric

#### 20-30 mm Hg Compression

Specify colors: BL-Black, BE-Beige, NU-Nude (no Nude in 3XL and 4XL) Open Toe in BL-Black and BE-Beige only Example: 9423-BL

#### 30-40 mm Hg Compression

Specify colors: BL-Black, BE-Beige Example: 9523-B

BL-Black and BE-Beige only Example: 9423-BL								
	Order #	Small	Medium	Large	X Large	2X Large	3X Large	4X Large
	20-30	9423	9425	9427	9428	9429	9430	9431
	Open*	OT - 9423	OT - 9425	OT - 9427	OT - 9428	OT - 9429	OT - 9430	OT - 9431
	30-40	9523	9525	9527	9528	9529	N/A	N/A
	Open*	OT - 9523	OT - 9525	OT - 9527	OT - 9528	OT - 9529	N/A	N/A
	Ankle	61/4" - 71/2"	7½" - 9"	9" - 10"	10" - 111/4"	111/4" - 121/2"	12½" - 13 ¾"	13¾" -15"
	Calf	10½" - 13"	11½" - 14"	12½" - 15"	13½" - 16"	14½" - 17"	16½" - 20½"	20" - 24"
	Length	11" - 15"	12" - 17"	13" - 18"	14" - 19"	16" - 21"	16" - 21"	16" - 21"

# Thigh High with Uni-Band

- Uni-Band Top for Men and Women with Silicone Beads.
- Available in Open or Closed Toe
- Traditional, Thicker, Non-Sheer Fabric

#### 20-30 mm Hg Compression

Specify colors: BL-Black, BE-Beige, NU-Nude Open Toe in BL-Black and BE-Beige only Example: 9433-BL

#### 30-40 mm Hg Compression

Specify colors: BL-Black, BE-Beige Example: 9533-BE

ĺ		Small	Medium	Large	X Large	2X Large
	20-30 Closed	9433	9435	9437	9438	9439
	20-30 Open*	OT - 9433	OT - 9435	OT - 9437	OT - 9438	OT - 9439
	30-40 Closed	9533	9535	9537	9538	9539
	30-40 Open*	OT - 9533	OT - 9535	OT - 9537	OT - 9538	OT - 9539
	Ankle	61/4" - 71/2"	7½" - 9"	8¾" - 10"	10" - 111/4"	11¼" - 12½"
	Calf	10½" - 13"	11½" - 14"	12½" - 15"	13½" - 16"	14½" - 17"
	Thigh	15" - 18"	17½" - 20½"	20" - 23"	22½" - 25½"	25" - 28"
	Length	23" - 26"	23" - 29"	25" - 31"	27" - 33"	35" - 41"



Open Toe

# **Pantyhose**

- Control Top
- · Reinforced Toe
- Traditional, Non-Sheer Fabric

#### 20-30 mm Hg Compression

Specify colors: BL - Black, BE - Beige

Example: 9443-BL

Order #	Small	Medium	Large	X Large	2X Large
Pantyhose	9443	9445	9447	9448	9449
Ankle	61/4" - 71/2"	7½" - 8¾"	8¾" - 10"	10" - 11¼"	11¼" - 12½"
Calf	10½" - 13"	11½" - 14"	12½" - 15"	13½" - 16"	14½" - 17"
Thigh	15" - 18"	17½" - 20½"	20" - 23"	22½" - 25½"	25" - 28"
Нір	Up to 40"	Up to 42"	Up to 44"	Up to 46"	Up to 52"
Waist	Up to 30"	Up to 331/2"	Up to 36"	Up to 38"	Up to 44"
Length	29" - 35"	31" - 37"	33" - 39"	35" - 41"	35" - 41"



## **Planograms**











#### Disclaimer

Some of our dealers bill medicare and private insurance as part of their business. Therefore Advanced Orthopaedics provides suggested HCPC codes for many products. However, PDAC periodically makes HCPCS coding changes. All changes and updates are posted on their website.

Therefore, we cannot guarantee coverage for a product under our suggested HCPCS codes and you are solely responsible for obtaining updated information from your carrier and properly coding and billing for Advanced Orthopaedics products and other products. Please visit our website for all current PDAC letters on file.

## **Product Return Policy**

Before returning any product, please call our Customer Service Department for return authorization. All returns must be accompanied with the authorization, invoice number and the date of purchase. No returns will be accepted after 30 days from the date of purchase. Please return all products in their original container.

All returns are subject to a restocking fee:

Return Policy
1 - 30 days
(after invoice date)

Restocking Fee 25%

1. This new Return Policy is effective as of July, 2021.

All credits are predicated on product inspection and approval for resell.
 Any product that does not meet these standards to be resold will not be credited.





#### **CUSTOMER SERVICE**

U.S. and International Headquarters 326 Main Street Harleysville, PA 19438 Toll Free: (800) 270-7074 Tel: (215) 256-9000

Fax: (215) 256-9006 orders@advanced-orthopaedics.com

#### **SALES & MARKETING**

Executive Offices
19 Shelter Cove Lane
Suite #207
Hilton Head Island, SC 29928
Tel: (843) 682-4400
Fax: (843) 682-4401

info@advanced-orthopaedics.com

#### **LOGISTICS MANAGEMENT**

Warehouse 5543 Columbia Blvd. Bloomsburg, PA 17815